


E faafefea ona aoina lau faata'ita'iga



1 Fufulu ou lima i le
vala'au.



2 Tatala le afifi o swab
ma 'alo'ese aua ne'i
pa'i i le pito o le
vavae.



3 Vilivili faatolu le pito
i luga o le swab
i puisu uma e lua.



4 Tuu le swab i totonu
o le fagu, faaū i lalo
le pito o le vavae.



5 Toe fufulu ou lima i se
vaila'au. Ua Mae'a!

E faafefea ona Lesitala le Pusa e fai ai Su'ega [Test Kit] o le COVID-19

1

Lesitala lau pusa-su'ega (poo le pusa-su'ega a lou alo) i luga o le initaneti i le color.com/covid/activate.

O le a mana'omia lau fa'aailoga faapitoa [barcode], lea e te maua i luga o lau fagu e aoina ai faata'ta'iga.



! **LAPATA'IGA:** E lē siakiina lau faata'ita'iga [sample] pe a lē lesitalaina lau pusa-su'ega.

2

Ina ia lesitalaina lau pusa-su'ega e mana'omia le faia o sau lesitala [account] ia Color.

1. Kiliki + Faapoopo tamaiti tausī
2. Faamolemole faatumu faamatalaga a le tamaitiiti o loo tausī.
3. Faatumu le numera ID ao'ga a le tamaititi tausī. Afai e te le silafia le ID a'oga, fesili i le aufaigaluega o le ofisa i luma i le a'oga.
4. O lenei laasaga e faamautinoa ai o loo matou feso'ota'iina lau lesitala [account] i le polokalame sa'o.

Mo Tamaiti a'oga e ta'i 14 pe i lalo ifo tausaga le matutua, e mana'omia e matua poo matua tausī ona lesitala se account e fa'aaogā ai le igoa o le matua/matua tausī ma tuātusi imeli, faatasi ai ma le ID o le Tamaititi a'oga e lesitala ai a latou pusa-su'ega.

Mo Tamaiti A'oga e ta'i 14 pe sili atu tausaga le matutua, e tatau ona fai ni a latou account e fa'aaogā ai o latou igoa ma tuātusi imeli, faatasi ai ma a latou ID a'oga e lesita ai a latou pusa-su'ega.

Mo le Aufaigaluega, e tatau ona fai ni a latou account e fa'aaogā ai a latou tuātusi imeli e feso'ota'i ma a latou a'oga.